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
Figure 16 2-lumen cystometry/pressure flow catheter

此導管用來測2壓力頻道Cystometry和Pressure-flow

- MMS 5702 2-lumen catheter (8 French). One for measuring the bladder pressure, the second for filling the bladder.
- MMS 5716 Paediatric 2-lumen catheter (6 French). One for measuring the bladder pressure, the second for filling the bladder.
- MMS 5415I Rectal PVC balloon catheter (4,5 French) with introducer for measuring the abdominal pressure.
- MMS 5349 2-lumen Cystometry catheter with conductance rings and modular telephone plug for Solar SHP/C 9 French

7.3 3-lumen Cystometry / Pressure-flow / UPP


Figure 17 3-lumen cystometry/pressure-flow/UPP catheter

此導管用來測3壓力頻道Cystometry, Pressure-flow, UPP和 Gynecology studies

MMS 5331 3-lumen catheter (9 French). One for measuring the bladder pressure, the second for filling the bladder and the third for measuring the urethra pressure.

MMS 5415I Rectal PVC balloon catheter (4,5 French) with introducer for measuring the abdominal pressure.

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MMS 302-07 3-lumen catheter for women (7 French), 20 cm.

MMS 302-09 3-lumen catheter for women (9 French), 20 cm.

MMS 312-07 3-lumen catheter for men (7 French), 40 cm.

MMS 312-09 3-lumen catheter for men (9 French), 40 cm.

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7.5 3-lumen UPP

- Put the DT-XX water transducers in the water transducer interface.
- Connect the transducer cables to the connectors on the interface.
- Place the water transducer holder at Symphysis Pubis level.
- Insert the bladder catheter and the rectal catheter.
- Connect the catheters to the water transducers.
- Flush the transducers/catheters to expel all air bubbles (see figures 21 and 22). Make sure all air is removed.
- Open the stop-cocks to atmosphere.
- Click the **Zero all** button to zero all pressure channels.
- Open the stop-cocks to the catheter.


Figure 21 Flush and zero with syringe

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WARNING Do not Slam and Gold Care no softly on pump. Do not switch on the roller pump if the C-clamp is removed. Never put fingers or any other body parts in the pump while it is running! Do not put any objects in the roller pump. The pump is designed to maintain a constant speed; the rolls of the pump will not be stopped easily.

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WARNING Silver pump. Do not switch on the roller pump if the pump cover is opened. Never put fingers or any other body parts in the pump while it is running! Do not put any objects in the roller pump. The pump is designed to maintain a constant speed; the rolls of the pump will not be stopped easily.

警示 確認保險絲(電源單元)的特別型式和數字以便更換 pump. The pump is designed to maintain a constant speed; the rolls of the pump will not be stopped easily.

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WARNING Do not hang weights > 3.5 kg on the hook of the Pump Infused volume Module. The infused volume transducer will be damaged.

警示 If the computer is connected to external volume devices (such as printer, modem, telephonelines,

WARNING Do not drop the flow meter. Do not apply excessive force on the flow meter (maximum 4 kg). This can permanently damage the flow meter combined with other equipment (for example Barostat or TMPD) make sure that the

requirements for safety. 1. Collateral standard requirements for safety) and/or IEC 60601-1 (Medical electrical equipment, Part 1: General requirements for safety).

WARNING Never attempt to open the remote control.

WARNING Do not attempt to open the remote control.

WARNING The digital cables connected to the Solar Main module (for example the flow meter) all have a connector with locking mechanism. Do not attempt disconnection by pulling the cable. Loosen the lock by pressing the lever on the top side of the connector.

WARNING Make sure you use a shielded cable to connect the printer.

WARNING Use only sterile filling liquids, disinfected or

WARNING Use only sterile filling liquids, disinfected or

should be seen as an aid. They are not substitutes for careful personal monitoring and an alert investigator.

WARNING Never turn off the Solar system if Windows is still active. You may loose data.

WARNING During infusion into the corpus cavernosum, extreme care must be taken that no air bubbles are infused into the patient's corpus cavernosum. Monitor the infused liquid constantly. Infusion of air can cause serious damage to the patient and is possibly lethal. If air bubbles are detected during the investigation, the pump must be stopped immediately.

WARNING High pressure can cause severe damage to the patient during cavernosometry. The pump will switch off on over-pressure only if the pressure is measured correctly. Instruct the patient to warn you if extreme pressure or pain is felt; then stop the pump immediately.

WARNING During a pre-test the pump will not stop on over-pressure (during cavernosometry)!

WARNING When using the catheter puller, do not clamp the catheter at the position of a pressure transducer, as this will damage the transducer permanently

NOTE Check the system regularly with the diagnostic software.

NOTE For cleaning of the catheters, we refer to the manufacturer's instructions.

WARNING Safety can be violated if, for example, the apparatus:

- shows visible damage.
- fails to perform the intended measurement.
- has been subjected to prolonged storage under unfavourable conditions.

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- has been subjected to severe stresses in transit.
- has been handled by unqualified personnel.

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Although precautions have been taken in this apparatus to prevent malfunctions due to abnormal environmental conditions and/or handling, such malfunctions can occur in some

T
y
p
e

B
F

E
q
u
i
p
m
e
n
t

Abnormal environmental conditions include, but are not limited to:

- any other condition as is specified in the manual.
- (electro-)magnetic fields caused by portable equipment (as telephones, computers), X-ray equipment, (radio-)transmitters, radar-equipment, etcetera.
- Electro Static Discharge.
- mains disturbances.
- extreme temperatures, pressures, humidity, dust, etcetera.

Abnormal handling includes, but is not limited to:

- any handling other than specified in the manual.
- rough handling.
- handling by unqualified staff.

Provided such or similar conditions do not occur, this apparatus will prove to be durable and reliable.

T
y
p
e

B

E
q
u
i
p
m
e
n
t

E
q
u
i
p
m
e
n
t

c
a
n

o
n
l
y

b
e

u
s
e
d

i
n
d
o
o
r

C
E

(
9
3
/

4
2
/
E
E
C
)
a
p
p
r
o
v
e
d

F
u
s
e

E
q
u
i
p
o
t
e
n
t
i
a
l
i
t
y

(
p
r

o
t
e
c
t
i
v
e

e
a
r
t
h
)

N
e
v
e
r

p
u
t

f
i
n
g
e
r
s

o
r

a
n
y

o
t
h

e
r

b
o
d
y

p
a
r
t
s

i
n

t
h
e

p
u
m
p

w
h
i
l
e

i
t

i
s

r
u
n
n
i
n
g
!

I
C
l
a
s
s

I
E
q
u
i
p
m
e
n
t

IPx1

P
r
o
t
e
c
t
e
d

a
g
a
i
n
s
t

d
r
i
p
p
i
n
g

w
a
t
e
r

A
t
t
e
n
t
i
o
n
,
c
o
n
s
u
l
t

t
h
e

m
a
n
u
a
l
s

Main module

USB ports
UP = Upstream from Solar main module to PC.
1..4 = Downstream from device to Solar main module.

Digital ports 1..8 I/O connections.

Pump Infused volume Module

Maximum weight infused volume module
3,5 kg.

Flush button pump filling tube.

Abort button for pump module (filling).
Only available on the Gold and the Gold Cavernosometry pump.

Direction of the pump module clockwise.
Only available on the Gold pumps.

Direction of the pump module counterclockwise.
Available on the Silver and Gold pumps.